Newborn Bathing and Nursery Care Model



Male









This training model is a life-size seamless manikin of a 0-4-week old infant, which has the weight and balance of the newborn as well as palpable anterior and posterior fontanel.

- 1. Life-like, life-size newborn model
- 2. Non-rigid neck flexibility
- 3. The fonticulus anterior and posterior are palpable.
- 4. The limbs can bend and stretch.
- 5. Waterproof: the model has a seamless, one-piece structure.
- 6. The umbilical cord is removable.

- Bathing
- Umbilical cord care
- Bottle-feeding
- Enema
- Ear and nose care
- Diaper and clothing changing
- · Physical measurement

- 1 newborn manikin
- 1 umbilical cord
- 1 feeding bottle
- 1 baby clothes 1 cloth cover

50 cm

Newborn Vital Signs Examination Simulator

M58



Product Supervision Keiko Takahashi, professor College of Nursing Science Fukui Prefecture University

Newborn Vital Signs Examination Simulator is an effective tool for nurses to learn how to take vital signs of a newborn.

- 1. Arterial pulse in wrists, neck and temples.
- 2. Rectal temperature can be checked.
- 3. Heart sounds are recorded from actual patients and synchronized with arterial pulses.
- 4. All vital signs can be controlled with the control unit.

- Vital Signs Examination
 - -Heart sounds
 - -Respiration
- -Pulses: superficial temporal, carotid and radial arteries (bilateral)
- -Body temperature (rectal thermometry)

- 1 newborn manikin
- 1 control box
- 1 thermometer 1 baby clothes

50 cm