# Difficult Airway Management Simulator -Demonstration Model-





# **Great addition as a demonstration tool for airway management**

Compact table-top design, anatomical accuracy and a variety of DAM settings. Ideal for demonstration of skills and devices.

### **Features**

- 1. Anatomically correct airway
- **2.** The incisors are removable when excessive force is applied.
- **3.** Successful tube tip placement can be confirmed by indicators.

## **Training Skills**

- Airway opening techniques (head tilt, jaw thrust)
- Bag-Valve-Mask ventilation
- Pre-intubation airway assessment
- "Sniffing position"
- Pressurization of external larynx to improve the laryngeal view
- Intraoral/Intranasal Intubation
- Use of oropharyngeal airway (OPA)
- Use of nasopharyngeal airway (NPA)
- Use of laryngeal mask airway
- Use of video laryngeal scope
- Confirmation of successful ventilation by indicators.
- Feedback of incorrect procedures including esophagus intubation and unilateral intubation
- Securing the tube in place with tape or Thomas™ endotracheal tube holder

### **Recommended Devices:**

- (Macintosh laryngoscope) blade size 3,4
- (Tracheal tube) internal diameter 7.0mm, 7.5mm
- (Laryngeal mask) Air-Q: size 3.5 LMA Supreme: size 4 I-qel: size 4

#### Specifications:

Specifications :	
Set Includes:	Manikin Size:
1 manikin	W45 x D20 x H25 cm
3 upper incisors	
1 lubricant	Replacement Parts:
1 syringe	(11393-010) 1 face mask
1 carrying bag	(11392-080) 2 replacement tongues II
1 instruction manual	(11392-010) 10 upper incisors



# Realistic trainings in insertion and management of a fiberscope


Anatomically correct trachea and bronchi provides realistic view through a scope, allowing recognition of bifurcations.

#### **Features**

1. Anatomically correct trachea and bronchi that include trachea bifurcation for segmental bronchi.

## **Training Skills**

- Insertion and management of a fiberscope
- Recognition of bifurcations by a fiberscope
- Airway opening techniques (head tilt, jaw thrust)
- "Sniffing position"



### Specifications:

Set Includes:	Manikin Size:
1 manikin	W45 x D20 x H25 cm
3 upper incisors	
1 lubricant	Replacement Parts:
1 syringe	(11393-010) 1 face mask
1 carrying bag	(11392-080) 2 replacement tongues II
1 instruction manual	(11392-010) 10 upper incisors