



Elastic external genitalia and realistic shape of vaginal canal realize a highly realistic training sessions on the flow of skills from the perineal protection to the delivery of placenta.

**Features**


1. Genital area made of innovative material that excels in elasticity, durability and resilience offers repetitive training opportunity.
2. Anatomically correct vaginal canal, ischial spine and pubic joint that facilitate training with anatomical understanding.
3. Elastic external genitalia, realistic shape of vaginal canal and the holding handle for the maternal body facilitate simulating a smooth delivery

**Training Skills**

- Delivery assistance for different delivery positions (dorsal, all-fours)
- Perineal protection
- Delivery of fetus (forceps delivery, vacuum extraction)
- Clamping, tying and cutting of the umbilical cord)
- Delivery of placenta
- Inspection of velamen

**Key Features**

**Realistic delivery of a baby with an elastic and resilient external genitalia**



**Skills training associated with anatomical understanding**

Position of the fetus against the pelvis can be observed under direct vision during the delivery.



**Two possible birth positions**

- Dorsal position
- All-fours position

**Specifications :**

<b>Set Includes:</b>		<b>Manikin Size:</b>	<b>Replacement Parts:</b>	
1 maternal body	1 talcum powder	W74 x D44 x H33.5 cm	(11415-030) 1 genitalia unit	(11416-010) 1 fetus
1 genitalia unit (for pelvic examination and delivery assistance)	5 velamen sheets		(11416-020) 1 placenta	(11416-030) 10 umbilical codes for omphalotomy
1 fetus	3 lubricant (150 ml)		(11416-040) 5 velamen sheets	(11415-020) 3 lubricant (150 ml)
1 placenta	1 vat		<b>Replacement Parts:</b>	
3 umbilical cords	1 instruction manual		(11417-020) 1 genitalia unit (perineorrhaphy) *skin included	