




Represents a perineal tear that occurs in a vaginal delivery. The innovative tissue material allows repetitive trainings.

## **Features**

- 1. Suture training for the first-degree perineal tear
- **2.** Genita area that excels in elasticity, durability and resilience for multiple training opportunity.
- **3.** Optional genitalia unit for training with the materbal body.

## **Training Skills**

- Perineorrhaphy (perineal repair)
  - ▶ Single interrupted suture
  - ▶ Vertical mattress suture



First-degree perineal tear



\* Suturing instruments are not included.



Genitalia skin can be mounted on the genitalia unit (optional part) for realistic training with the maternal body.

## Specifications:

| Set Includes:        | Manikin Size:        | Replacement Parts:                            |
|----------------------|----------------------|---|
| 1 ganitalia skin     | W16 x D23 x H20.5 cm | (11417-010) 2 genitalia skin (perineorrhaphy) |
| 1 supporting holder  |                      | Optional Parts:                               |
| 1 talcum powder      |                      | (11417-020) Genitalia unit (frame, skin)      |
| 1 instruction manual |                      |   |