Peripheral Venous Catheter Placement Simulator

MW9



Peripheral Venous Catheter Placement Simulator provides full procedural skills from injection to fixing the IV cannula.

Features

- 1. All training skills for the peripheral venous route management
- 2. Two puncture sites: median antebrachial vein and dorsal vein
- 3. Practice for the extension of hand skin and the angiopressure management
- 4. Durable puncture pad with IV cannula

Training Skills

- Wrapping with tourniquet
- Confirmation of puncture site
- Sterilization of puncture site
- Puncture with IV cannula
- · Confirmation of back flow in puncture
- Angiopressure management and decannulation
- Setting of influsion tube
- Confirmation of natural instillation
- Fixing of the puncture site
- Injection of medical solution from injection sub port







Set Includes

- 1 arm model
- 2 injection pads for dorsal vein of hand
- 2 injection pads for median antebrachial
- 1 bottle for simulated blood
- 1 stand for bottle of simulated blood
- 1 infusion set
- 1 infusion bag
- 1 syringe (50ml)
- 1 plastic beaker
- 10 simulated blood sticks (swab type)
- 1 instruction manual

Manikin Size • • • Replacement Parts

- W52 x D10 x H10 cm
- 11388-200
- 2 injection pads for median antebrachial vein
- 11388-300
- 2 injection pads for dorsal vein of hand
- 11388-400
- 10 simulated blood sticks (swab type)