

M86

Child Patient Care Model

5-6 years old

Instruction Manual

Index

This training model is a life-size seamless manikin of a 5-6-year old child.

The simulator is designed for training in nursing care.

Do not use for improper purpose.

Please read the instruction carefully before use.

Specifications

1. Life-like, life-size 5-6 year old model
2. The model is designed to provide the feeling of approx. 15kg weight for patient handling practice.
(15 kg is not the actual weight)

Set includes

1 manikin

1 set of pajamas

DOs and DON'Ts

DOs

Handle with care

The model is made of special composition of soft resin. Please handle it with the utmost care like an actual baby at all times. Store the manikin with its clothes on, covered by the supplied cloth, away from heat, moisture and direct sunlight.

DON'Ts

The model is not designed for bathing. **NEVER immerse the manikin in water.**

Don't mark on the model with pen or leave any printed materials in contacted with their surface.

Ink marks on the models will not be removable.

KYOTO KAGAKU CONTACT & ORDERING INFORMATION

Inquiries & Ordering:

Telephone +81-75-605-2510

Facsimile: +81-75-605-2519

URL: <http://www.kyotokagaku.com/>

Email: rw-kyoto@kyotokagaku.co.jp

Head Office:

Kyoto Kagaku Co. Ltd.

15 Kitanechoya-cho, Fushimi-ku

Kyoto, Japan 612-8388

USA & CANADA:

Kyoto Kagaku Co. Ltd. LA Office

13109 Lomita Boulevard, Torrance, CA90505, USA

Telephone 310-325-8860

Facsimile: 310-325-8867

MEM Training Centre:

Kyotokagaku Co., Ltd. Tokyo Branch

5-20-4 Koishikawa Bunkyo-ku Tokyo, Japan 112-0002